


Knee Revision Set

Used for total knee revision



PRODUCT NO'S:	
MS5470-00 [Complete Set]	
Individual Instruments: MS5470-08 [8mm Chisel]	MS5472-08 [8mm Offset Cement Removal Chisel] Osteotome Dimensions: 8mm Wide x 12mm Long Blade Length: 2.375" Overall Length: 7.375"
Osteotome Width: 8mm Blade Length: 2.375" Overall Length: 7.375"	MS 5474-06 [6mm Notched Cement Removal Chisel] Osteotome Width: 6mm Blade Length: 2.625" Overall Length: 7.375"
MS5470-11 [11mm Chisel] Osteotome Width: 11mm Blade Length: 2.375" Overall Length: 7.375"	MS5475-08 [8mm Implant Remover] Diameter: 8mm Blade Length: 2.625" Overall Length: 7.375"
MS5470-20 [20mm Chisel] Osteotome Width: 20mm Blade Length: 2.375" Overall Length: 7.375"	MS 5470-CASE [Case Only]

MED SURGICAL INSTRUMENTS 724, Thornhill Drive Coliegeville. PA 19426.

TEL: +1 610-983-3073 | Email: info@medsurgicalusa.com | Web: www.medsurgicalusa.com

